RESMED

CLINICIAN'S GUIDE

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MODE D'EMPLOI

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BRUKSANVISNING

RESMED

USER'S GUIDE (ENGLISH)	1
GEBRAUCHSANWEISUNG (DEUTSCH)	13
MODE D'EMPLOI (FRANÇAIS)	24
GUIDA PER L'UTENTE (ITALIANO)	35
GUÍA DEL USUARIO (ESPAÑOL)	46
GUIA DO USUÁRIO (PORTUGUÊS)	57
BRUKSANVISNING (SVENSKA)	68

CONTRAINDICATIONS	2
WARNINGS	2
PRECAUTIONS	2
ADDITIONAL INFORMATION	3
PARTS OF THE MASK	4
FITTING THE MASK	4
REMOVING THE MASK	6
QUICK RELEASE TAB	6
CLEANING THE MASK IN THE HOME	6
CLEANING AND DISINFECTING/STERILISING THE MASK BETWEEN PATIENTS	7
REASSEMBLING THE MASK	9
MASK CHARACTERISTICS	10
ACCESSORIES AND REPLACEMENT PARTS	11

MIRAGE™ NV FULL FACE MASK SERIES 2

The MIRAGE™ NV FULL FACE MASK SERIES 2 system (non-vented) is only indicated for the treatment of adult patients using non-invasive ventilator systems that have an active exhaust.

CONTRAINDICATIONS

- The MIRAGE NV FULL FACE MASK SERIES 2 must not be used without appropriately qualified supervision either when the patient is at risk of not responding to severe hypoxia or hypercapnia by arousal, or when the patient, having aroused, is unable to remove the mask.

WARNINGS

- The mask must not be used on patients with impaired laryngeal reflexes or other conditions predisposing to aspiration in the event of regurgitation or vomiting.
- The mask should not be worn unless the ventilator system is turned on and operating properly.

- At a fixed rate of supplemental oxygen flow (if used), the inhaled oxygen concentration will vary depending on the pressure settings, patient breathing pattern, mask selection and the leak rate.
- The mask must not be used if the patient is known to be allergic to any of the following mask component materials: polycarbonate, silicone rubber, polyurethane rubber, nylon, polyester and cotton.

PRECAUTIONS

- The mask is not for use with ventilators which have a combination of:
 - a) little or no ability to compensate for mask leak and variation in mask leak (eg volume-cycled ventilators with no compensation for leak, or pressure-cycled ventilators which are unable to deliver adequate extra flow in the presence of increased leak); and
 - b) no expired ventilation alarms.

If a ventilator with poor leak compensation is used, special reliance is placed on the expired ventilation alarm.

- As with all mask ventilation systems, significant leak may occur between the mask and the patient's face. The leak may vary depending on such factors as the patient's position and sleep stage. Variation in leak may significantly affect patient ventilation, potentially producing profound hypoventilation or hyperventilation. The magnitude of this effect depends on the ventilator's ability to compensate for leak.
- Leak or variation in the leak may cause mistriggering of the ventilator.
- Increased face-to-mask leak will reduce effective deadspace by washing out expired gases. In a pressure-cycled ventilator which tolerates moderate leak, the effect may be a reduction in arterial PCO_2 with no change in ventilator settings.
- In determining the initial ventilator settings, and in particular the tidal volume for volume-cycled ventilators, make appropriate allowance for mask dead space and leak. Refer to the Mask Characteristics section of these instructions for more information.
- Refer to the relevant ventilator system operating instructions for detailed information concerning connecting the mask.

ADDITIONAL INFORMATION

- ResMed has tinted the **MIRAGE NV FULL FACE MASK SERIES 2** blue to help you identify it as a non-vented mask for use with active-exhaust non-invasive ventilator systems only.
- The mask has a standard female 22mm tapered fitting (ISO5356-1) for external connection to the ventilator.

PARTS OF THE MASK


Figure 1: Mask Parts


Figure 2: Headgear

FITTING THE MASK

The MIRAGE NV FULL FACE MASK SERIES 2 is supplied fully assembled.

Fitting the mask correctly is essential to optimum therapy delivery. Use the ResMed Mask Fitting Template to help you choose the best mask size for your patient.

- 1 Connect the mask to the ventilator system by following the ventilator system operating instructions and/or physician's prescription.
- 2 Seat the mask cushion over the patient's nose and mouth. On most patients, the bottom of the cushion should rest below the lower lip and the top of the cushion on the bridge of the nose, between the eyes.
- 3 Pull the headgear over the patient's head. Make sure the side lower Velcro® straps pass below the patient's ears, and that the side upper Velcro straps pass over their forehead.

- 4 Tighten the side lower (Fig. 3) and then the side upper Velcro straps (Fig. 4) to achieve a secure and comfortable fit.


Figure: 3


Figure: 4

- 5 Adjust the top Velcro strap (Fig. 5) with very slight tension to add stability to the fit of the mask. The mask should be secure but not overly tight.


Figure: 5

- 6 The fitted mask should look like this (Fig. 6):


Figure: 6

Note: The top of the mask frame should not touch the patient's forehead. In most cases, the top of the mask frame should be some distance away from the patient's forehead.

- 7 Turn the ventilator on. If you detect air leaks around the cushion, reposition the mask on the patient's face to achieve a better seal. You may also need to readjust the Velcro straps.

Note: When the ventilator system is delivering positive pressure, the Mirage cushion inflates to float gently on the patient's face, forming a secure and comfortable seal without the need for tight headgear. Pulling the headgear too tightly may reduce the effectiveness of the mask.

REMOVING THE MASK

To remove the MIRAGE NV FULL FACE MASK SERIES 2, release the lower Velcro straps and pull the mask and headgear up over the patient's head.

Note: The Quick Release Tab should be used to remove the mask only in the event of an emergency.

QUICK RELEASE TAB

The Quick Release Tab (Fig. 7) allows the mask to come away from the face quickly in the event of an emergency.

Simply pull the Quick Release Tab and the rear Velcro opening of the headgear separates, allowing you to remove the mask from the face.


Figure: 7

Note: The Quick Release Tab should be used to remove the mask only in the event of an emergency. Overuse of the Quick Release Tab will result in the deterioration of the Velcro opening.

CLEANING THE MASK IN THE HOME

Note: Wash the headgear before first time use, as the dye may run. It is important to always handwash the headgear.

- 1 Detach the mask from the ventilator and separate the headgear, luer lock port caps, cushion clip and mask cushion.
- 2 Handwash the headgear in warm, soapy water. Use a mild detergent such as pure soap. Rinse well and allow it to air dry out of direct sunlight. Do this at least once per week.
- 3 Handwash the mask components in warm, soapy water. Use a mild detergent such as pure soap. Rinse the components well and allow them to air dry out of direct sunlight. Do this after each use.

PRECAUTION

- Do not use bleach-, chlorine-, alcohol-, or aromatic-based solutions, moisturising or antibacterial soaps, or scented oils to clean any part of the mask system, including headgear. These solutions may cause damage and reduce the life of the product.
- Do not expose any part of the mask system or tubing to direct sunlight as it may cause the product to deteriorate.

CLEANING AND DISINFECTING/STERILISING THE MASK BETWEEN PATIENTS

The following are the recommended and validated procedures for cleaning, disinfection and sterilisation of the MIRAGE NV FULL FACE MASK SERIES 2. We advise that each healthcare facility consult its own procedures before carrying out the following instructions.

WARNING : ResMed cannot give any assurance that deviations from the procedures listed below, and their effect on the performance of the product, will be acceptable.

- 1 Detach the mask from the ventilator. Disassemble the mask by separating the headgear, luer lock port caps, cushion clip and mask cushion.
- 2 After disassembly, thoroughly clean the mask components according to the following cleaning procedure:

CLEANING PRIOR TO DISINFECTION OR STERILISATION

- a Handwash the headgear in warm water using a mild detergent, such as pure soap. Rinse well and allow it to air dry out of direct sunlight.

- b Clean the individual parts of the mask with a soft bristle brush for one minute while soaking them in detergent. Pay particular attention to all crevices and cavities.
- c Rinse the components twice by agitating them vigorously in one litre of water and allow them to air dry out of direct sunlight.

PRECAUTION

- Do not expose any part of the mask system or tubing to direct sunlight as it may cause the product to deteriorate.
 - Do not disinfect or sterilise the headgear. This part cannot be disinfected or sterilised.
- 3 After disassembly and cleaning, disinfect the reusable mask components according to the following procedure for chemical or thermal disinfection.

OR

After disassembly and cleaning, sterilise the mask according to the following procedure for the STERRAD® Sterilisation Process.

Note: The mask must be disassembled and cleaned prior to disinfection or sterilisation.

OPTION 1: CHEMICAL DISINFECTION

- a** Soak the reusable mask components in a commercially available solution of glutaraldehyde 3.4% (eg CIDEX Plus®), peracetic acid 0.08% (eg CIDEX PA®) or orthophthalaldehyde 0.55% (eg CIDEX OPA®) as described in the manufacturer's instructions.
- b** Rinse the mask components as described in the manufacturer's instructions and allow them to air dry out of direct sunlight.
- c** Perform a visual inspection of each mask component. If any visible deterioration of a mask component is apparent (eg cracking, tears etc), the mask component should be discarded and replaced. Yellow discolouration of the cushion may occur. This discolouration is acceptable.

Testing by ResMed has shown that the mask components can withstand 15 cycles using solutions of glutaraldehyde 3.4% (eg CIDEX Plus®), peracetic acid 0.08% (eg CIDEX PA®) or orthophthalaldehyde 0.55% (eg CIDEX OPA®) when following the manufacturer's instructions.

OPTION 2: THERMAL DISINFECTION

- a** Using a certified hot water disinfection system, soak the mask components for one of the following times:
 - 70°C for 100 minutes or
 - 75°C for 30 minutes or
 - 80°C for 10 minutes.
 - 90°C for 1 minute.
- b** Perform a visual inspection of each mask component. If any visible deterioration of a mask component is apparent (eg cracking, tears etc), the mask component should be discarded and replaced. Yellow discolouration of the cushion may occur. This discolouration is acceptable.

Testing has shown that the mask components can withstand 15 cycles of disinfection using the above time-temperature combinations.

OPTION 3: STERRAD® 100S STERILISATION PROCESS

- a** Sterilise the mask components as described in the manufacturer's instructions for the STERRAD 100S Sterilisation System.

- b** Perform a visual inspection of each mask component. If any visible deterioration of a mask component is apparent (eg cracking, tears etc), the mask component should be discarded and replaced. Yellow discolouration of the cushion may occur. This discolouration is acceptable.

Testing by ResMed has shown that the mask components can withstand 15 cycles of the STERRAD 100S Sterilisation System.

WARNING: MIRAGE NV FULL FACE MASK SERIES 2 components should not be subjected to autoclave or ethylene-oxide gas sterilisation.

- 4** Reassemble the mask as described in the “Reassembling the Mask” section on page 9.

Note: If a healthcare facility requires an additional disinfection or sterilisation process after reassembly, the number of validated cycles is halved.

REASSEMBLING THE MASK

Ensure that all parts of the mask system, including headgear and Quick Release Tab mechanism, are clean, dry and free from damage before reassembling. Replace any affected part if it shows signs of wear or damage.

- 1** Twist the luer lock port caps onto the mask ports (Fig. 8).


Figure: 8

- 2** Fit the grooved edge of the mask cushion to the mask frame. Fit the top section first, then work your way around the mask, pushing the cushion onto the frame, ensuring a good seal has been achieved.

Then push the cushion clip over the cushion onto the mask frame, ensuring it clicks into place in all locations (Fig. 9).


Figure: 9

- 3 Thread the straps of the headgear through the slots on the mask and fasten the Velcro. The headgear label should face outwards and be near your neck when fitted.

MASK CHARACTERISTICS

MASK DEAD SPACE

Small mask	145 mL
Medium mask	200 mL
Large mask	230 mL

Note: These values were obtained using a facial profile and measuring the volume from the surface of the face to the end of the tapered fitting. The data is provided as an indication only. Actual volumes may vary depending on individual facial features.

MASK IMPEDANCE DATA

At a flow of 60 L/min the difference in pressure between the end of the tapered fitting and inside the mask is less than 0.5 cmH₂O.

MASK LEAK

The leak from the mask, excluding leak between the face and the mask, is less than 2.0 L/min at 30 cmH₂O (maximum pressure).

ACCESSORIES AND REPLACEMENT PARTS

Description	Size	Quantity	Part
		per pack	Number
Cushion	Small, standard	1	16604
Cushion	Small, shallow	1	16671
Cushion	Medium, standard	1	16605
Cushion	Medium, shallow	1	16672
Cushion	Large, standard	1	16606
Cushion	Large, shallow	1	16673
Cushion clip	Small	1	16674
Cushion clip	Medium	1	16675
Cushion clip	Large	1	16676
Headgear		1	16716
Luer lock port caps		2	16681
Luer lock port caps		10	16682
Mask frame	Small	1	16713
Mask frame	Medium	1	16714
Mask frame	Large	1	16715
Mask fitting template		1	16638

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© 2002 ResMed Ltd.
Velcro is a registered trademark of Velcro Industries B.V. CIDEX and STERRAD are registered trademarks of Johnson & Johnson
Mirage is a trademark of ResMed Ltd.